

High Capacity Ring Mill




- Maximum sample of 1500g
- Output particle size -75microns
- Air lifting assistance
- Full safety interlocking system
- Operate below 80db
- Carbon or Chrome Steel Head available



RM2000

The ROCKLABS RM2000 is designed to process large samples with a maximum sample weight of 1500g and a final particle size of -75microns. Fully constructed with metal for long lasting and powder coated finish for ease of cleaning.

The RM2000 has a pneumatically operated lifting arm for lifting and lowering the Head during loading and unloading of each process.



SOCACHIM



SPECIALIST IN SAMPLE PREPARATION

www.socachim.com

SOCACHIM, 282 Av. de Roodebeek, Brussels, Belgium.
Phone: +32 2 762 77 12. Fax: +32 2 762 55 07. Email: info@socachim.com

